







## Page 4 of 4 MEDICAL DEVICE WITH RADIOPACITY Verivada Chandrasekaran et al. 10527-410002

Sputter Clean	ning Parameters	Deposition S	Deposition Sputtering Parameters	Deposition	Deposition Parameters
(Ion	Beam)	(I)	(Ion Beam)	(Electron Beam Evaporation)	Evaporation)
Gas type	Argon	Gas type	Argon	Pressure	$\leq 1 \times 10^{-4} \text{ Torr}$
Flow	2.5-3.5 sccm	Flow	2.5-3.5 sccm	Deposition rate	3.5-3.5 Å/sec
Energy	1-2 keV	Energy	1-3 keV	Energy	8 keV
Power	58-65 Watts	Power	58-65 Watts	Power	2.9 kW
Current	10-35 mA	Current	15-35 mA	Current	0.375 A
Time	20-30 min				

## FIG. 7

,					
Sputter Clean	ning Parameters	Deposition S	Deposition Sputtering Parameters	Deposition	Deposition Parameters
(Ion Be	Beam)	D	Ion Beam)	(Electron Beam Evaporation)	n Evaporation)
Gas type	Argon	Gas type	Argon	Pressure	3.5x10 <sup>-5</sup> Torr
Flow	2 sccm	Flow	2 sccm	Deposition rate	4-5 Å/sec
Energy	2 keV	Energy	2 keV	Energy	8 keV
Power	60 Watts	Power	80 Watts	Power	1.76 kW
Current	15 mA	Current	18 mA	Current	0.226 A
Time .	20 min				

## FIG. 8

Sputter Cleaning	ning Parameters	Deposition S	Deposition Sputtering Parameters	Deposition	Deposition Parameters
el (Ion Be	Beam)		(Ion Beam)	(Electron Bear	(Electron Beam Evaporation)
Gas type	Argon	Gas type	Argon	Pressure	1.2x10 <sup>-4</sup> Torr
Flow	3 sccm	Flow	2 sccm	Deposition rate 2-3 Å/sec	2-3 Å/sec
Energy	1 keV	Energy	1 keV	Energy	7.5 keV
Power	60 Watts	Power	60 Watts	Power	2.0 kW
Current	25 mA	Current	32 mA	Current	0.15 A
Time	30 min				

FIG 9